



TECHNICAL DATA SHEET(TDS)
Alizarin Red Indicator Solution
Product Specifications (CAS# 130-22-3)

Test	Specification	Typical Results
pH Transition Interval	pH 4.6 (Yellow) to pH 6.0 (Pink)	Pass

Appropriate/legal use of this product is the responsibility of the user.

[Alizarin Red Indicator Solution](#)

Phone: 512-668-9918, Fax: 512-886-4008, E-mail: customercare@laballey.com
www.laballey.com

12501 Pauls Valley Road, Suite A, Austin, Texas 78737
© 2026 Lab Alley